

VistA

Audiometric Exam Module User Manual

Patch ACKQ*3.0*3

June 2003

Department of Veterans Affairs
VistA Health Systems Design & Development

Preface

Purpose

ACKQ *3.0*3 was designed to use a consistent, event-driven windows style clinical user interface to provide Audiologists and their staff with an easy way to enter, display, store and utilize the information obtained during the Audiometric exam of a patient.

Audience

The information in this manual is intended to aid practitioners in Audiology and Speech Pathology Service (ASPS) in the use of this software.

Benefits

Prior to this patch, there was no VA-wide availability of a way for practitioners to enter, store and view audiograms. Nor was there a central database for those readings. This patch introduces tools for these data entry, storage (both locally and centrally), and display capabilities. Upon structured entry of clinical hearing loss data through this package, a practitioner can immediately view an audiogram display and print a standard VA form 10-2364. An audiogram can also be printed for hard copy records or copied to the Windows clipboard for inclusion in other electronic documents.

Completed and signed audiograms are stored within the local *Vista* system. They are also transmitted electronically to the Denver Distribution Center (DDC), where they are stored centrally and are available to facilitate the ordering and maintenance of hearing aids and accessories. This audiometric information is subsequently available whenever a practitioner places a REMOTE ORDER ENTRY SYSTEM (ROES) order.

Related Manuals

Audiometric Module ACKQ*3.0*3 Implementation Guide
Audiometric Module ACKQ*3.0*3 Technical Manual
Audiometric Module ACKQ*3.0*3 Security Manual

This page left intentionally blank

Table Of Contents

Preface.....	i
Purpose.....	i
Audience.....	i
Benefits.....	i
Related Manuals.....	i
Table Of Contents.....	iii
Introduction.....	1
1. General Conventions Used.....	1
2. Applications.....	1
Audiometric Exam Enter/Edit.....	2
1. Enter/Edit Conventions.....	2
2. Accessing the Enter/Edit.....	3
Invoking From CPRS:.....	4
Invoking From the Desktop:.....	4
Continuing On:.....	5
3. Audiometry Entry Tab.....	7
4. Pure Tones Tab.....	10
5. Speech Audiometry Tab.....	12
Acoustic Immittance Tab.....	14
The Audiogram Display.....	15
1. Accessing the Audiogram Display.....	15
2. Display Conventions.....	19
3. Main View.....	21
Audiogram in overlapping view.....	21
4. Separate View.....	23
5. Viewing 10-2364.....	25
Glossary.....	27
Acronyms.....	27
Appendices.....	29
A: Determination of Series Values for Display.....	29
B: Calculation Of PB MAX And PI/PB.....	30
C: If You Have Access To Multiple Broker Environments.....	30
D: Calculation of Pure Tone Averages.....	31
E: VA FileMan Date/time Formats.....	31
Index.....	33

This page left intentionally blank

Introduction

1. General Conventions Used

All examples shown in this document are of fictitious patients and information and are meant only to show display features, not actual readings.

Functionality included in ACKQ*3.0*3 is implemented through two software applications, one providing primarily data entry tools and one providing primarily display tools. These applications use a Microsoft Windows-style graphical interface with typical mouse navigation and field selection. End users can also access editable fields, tab pages and menu options through the use of short-cut key combinations (using the ALT key in combination with some other key). The key to use is shown underlined and capitalized (i.e. Print). Occasionally the user may need to hold down the ALT key to activate the appropriate underlining of characters.

Additional conventions specific to each application are described in the corresponding sections of this document.

2. Applications

AUDIOGRAM DISPLAY [ACKQROES3]

AUDIOMETRIC EXAM ENTER/EDIT [ACKQROES3E]

Audiometric Exam Enter/Edit






1. Enter/Edit Conventions


The material on the following pages describes procedures and functional features applying to the navigation and data entry of a patient's audiometric readings. Below are related conventions specific to this application.


The application window includes multiple tab pages. Be sure that the window is sized large enough to see the tabs at the top and the buttons at the bottom of the pages.


All fields that are disabled (cannot be edited) have a grayed appearance. All tab stops have been arranged so the user can proceed in a logical sequence through the editable fields by using the Tab key. When duplicate data entry fields exist for both the right and left ears, keyboard navigation generally progresses through all of the right ear fields first, followed by all of the left ear fields. Masking level fields are editable only if the threshold value entered for the corresponding frequency includes the masking symbol (*). Without that symbol, Tab-key navigation proceeds across the page from threshold to threshold, skipping the masking level. Indicate 'No Response' by entering the plus (+) symbol.


At any time the user can bypass the Tab-key sequence by using the mouse to click in the desired field. Pausing the cursor over an editable field displays a hint with the acceptable range of values. Inactivated or disabled fields are not reachable with either the Tab key or the mouse.


The , , , , and  buttons appear on all tab pages.

The  button is disabled until an editable field has been changed. When this button is selected, changes made within the currently viewable tab page are stored.

The  button prints the currently viewable page.

If the  button is clicked from any tab, the audiogram display for the selected record is invoked. It is not necessary for the record to be signed in order to display the audiogram. See the [Audiogram Display](#) section for further information.


The  button displays help text specific to the currently viewable page.

The  button exits from the *Audiometric Exam Enter/Edit* application.

2. Accessing the Enter/Edit

Typically, this application can be invoked from either CPRS or a stand-alone desktop shortcut (see instructions below). If either of these methods is desired but is not available, contact your facility's IRM Service to have the necessary installation or setup procedures completed.

When invoking the application, if your facility does not have single signon enabled, you may need to log in with your local *VistA* Access and Verify codes. A form similar to the following will appear.




The image shows a Windows-style dialog box titled "VISTA Sign-on". The main text area contains a message surrounded by asterisks: "You are online with the Denver Distribution Center 'DEVELOPMENT' system. This system is available as a non-production simulated VistA environment for application development." Below the text is a stylized logo consisting of interlocking geometric shapes. To the right of the logo are two input fields: "Access Code:" followed by a field containing six asterisks, and "Verify Code:" followed by a field containing eight asterisks. To the right of these fields are two buttons: "OK" with a green checkmark icon and "Cancel" with a red X icon. At the bottom of the dialog, there is a status bar with four fields: "Server: AMMDDCDHC2", "Volume: TNR", "UCI: DEV", and "Port: //./nul:411".

In some cases, an end user may be set up for access to multiple broker environments. If that is the case, refer to [Appendix C](#) for further instructions.

Invoking From CPRS:

From the *Cover Sheet* in the CPRS application, click on the **T**ools menu. *Audiogram Edit* should be one of the options on that menu. Since a patient has already been selected in CPRS, you will not have to choose one.


At this point you will see the following question:



These instructions continue in the section [Continuing On](#).


Invoking From the Desktop:

Double-click the desktop icon for the *Audiogram Edit* application. After logging in (if necessary). At this point you will see the following question:




Your response will indicate whether the lookup for the patient is done in the PATIENT file or from the AUDIOMETRIC EXAM file.

If an existing audiogram is being entered, an existing exam will be selected from the AUDIOMETRIC EXAM file. Enter a few characters of the patients name in the following box.




Continuing On:

If you are not entering a new record and more than one record exists for the patient, you will need to select the specific exam from a screen similar to the following, and continue the instructions in the [Entry Tab Section](#).




If you are entering a new record you will see the following prompt:



Note: Exam DATE *AND* TIME must be included in this space for a new record to be created.

You may enter any of the standard VA FileMan date/time formats (see [Appendix E](#)). At this point, a new record is set up in the file and the age of the patient is calculated.


You will then see the following verification message:



When this box is closed the introductory tab page is shown, and you will continue on to the *Audiometry Entry* tab, as shown below.

3. Audiometry Entry Tab

The following message may appear before the application window is opened to the initial tab page. If so, click *OK* to continue, and adjust as needed to see all parts of the form.



Continue to the next page.

This page collects and displays basic patient and visit information. The DATE/TIME OF VISIT, PATIENT, AGE AT VISIT, and the VA ELIGIBILITY STATUS are calculated when the record is created and are not editable.

Audiometric Data for ZZTESTPATIENT,PAMS(0140) as seen on JUN 18, 2003

Audiogram Entry | Pure Tones | Speech Audiometry | Acoustic Immittance

DATE /TIME OF VISIT: JUN 18, 2003@08:40:10

PATIENT: ZZTESTPATIENT,PAMS

EXAMINING AUDIOLOGIST: URRUTIA,PAM

REFERRAL SOURCE: AUDIOLOGY

AGE AT VISIT: 63

VA ELIGIBILITY STATUS: Not Found

TYPE OF VISIT:

- ☐ C&P
- ☐ Audiological Assessment
- ☒ Other (check, then right click to edit)

(Entering a date here LOCKS the entire record and sends it to the DBC)


DATE SIGNED: [] Delete Record

Save Changes View Audiogram Print Help <?> eXit Program

The EXAMINING AUDIOLOGIST field is selectable from the New Person file (#200). Both the name and the title of the person selected are displayed on the audiogram. Entry of the first few characters of the person's last name will assist in the lookup.


The REFERRAL SOURCE is selectable from the Hospital Location file (#44). Again, putting in the first few letters of the location will assist in the lookup. Naming for these locations may vary from facility to facility.

If "Other" is selected from the TYPE OF VISIT radio buttons, a free text entry of 2 to 25 characters is required and stored with the record.

A screenshot of a Windows-style dialog box. The title bar is blue and contains the text "Enter 2 to 25 characters" followed by a close button (X). The main area of the dialog is light gray. It contains a label "Type of Visit" in a standard font. Below the label is a white text input field with the placeholder text "sample type of visit". At the bottom of the dialog, there are two buttons: "OK" and "Cancel", both in a standard font.

The VA ELIGIBILITY STATUS is automatically taken from the Patient file(#2) and is: "Verified", "Pending" or "Pending Re-Verification".

The record may be deleted as long as the DATE SIGNED field is empty. DATE SIGNED refers to the date that the information is verified as accurate, and once a date is entered into that field, the record is considered final and approved and the values are transmitted to the DDC database.

All fields become unavailable for further editing and the  button disappears.

4. Pure Tones Tab

Click on the *Pure Tones* tab or enter ALT-T to make the *Pure Tones* tab active. This tab allows for entry of pure tone threshold values across a standard range of frequencies as a measure of hearing loss.

Audiometric Data for ZZTESTPATIENT,PAMS(0140) as seen on JUN 18, 2003

Pure Tones

AIR CONDUCTION

	Hz	250	500	750	1000	1500	2000	3000	4000	6000	8000
RIGHT											
Initial											
Retest/Masked			60*+		70*	60	105*+	60*	69*	105+	100*+
Masking Level											
Initial			70		55		110				110
Final					65						
LEFT											
Initial											
Retest/Masked		20	70	40	50*+	60	60	70	70		90*+
Masking Level											
Initial					55						100
Final											

BONE CONDUCTION

	Hz	250	500	750	1000	2000	3000	4000
RIGHT								
Initial			30					
Retest/Masked			40*+		40	50	60*	60*+
Masking Level								
Initial			45					65
Final								
LEFT								
Initial			40	50	60	70	40	
Retest/Masked								45*+
Masking Level								
Initial			40					50
Final			40					

INTER-TEST CONSISTENCY

RIGHT: GOOD (dropdown) LEFT: POOR (dropdown)

PURE TONE AVERAGES

	RIGHT	LEFT
Two	65	55
Three	78	60
Four	76	62

SRT (Effective Masking Level)

	RIGHT	LEFT
Initial	60	70
Retest	65	70*
(Initial)		75
(Final)		

Buttons: Save Changes, View Audiogram, Print, Help <?>, eXit Program

Threshold values for each frequency tested should be entered into the corresponding input field. If the cursor is paused over an editable field, a hint with acceptable values is displayed. Along with numerical values, the symbols '*' for 'masked' and '+' for 'no response' can be entered in the Retest/Masked. Disabled fields and calculated fields that do not require data entry are displayed in gray.

Masking level fields are editable only if the repeat threshold value entered for the corresponding frequency includes the masking symbol (*). Without that symbol, Tab-key navigation proceeds across the page to each successive frequency, skipping the masking level field.

[Pure Tone Averages](#) (PTA) are automatically calculated when sufficient information is entered. See Appendix D for further explanation.

5. Speech Audiometry Tab

Click on the Speech Audiometry tab or enter ALT-A to make the Speech Audiometry tab active. This tab displays the Pure Tone Averages (see [Appendix D](#)) from the Pure Tone tab and allows for entry of comfort levels, effective masking levels and word recognition information. The PB Max and PI/PB fields are calculated when sufficient information is entered. See [Appendix B](#) for further information.

Audiometric Data for ZZTESTPATIENT,PAMS(0140) as seen on JUN 18, 2003

SPEECH COMFORT LEVELS

	MCL	UCL
RIGHT		
LEFT		

PURE TONE AVERAGES

	Two	Three	Four
RIGHT	65	78	76
LEFT	55	60	62

SRT

	Initial	Retest
Right	60	65
Left	70	70*

EFFECTIVE MASKING LEVELS

	Initial	Final	PB Max	RI
Right			80	0.06
Left	75		80	

WORD RECOGNITION TESTING

RIGHT EAR WORD RECOGNITION

	1	2	3	4	5
%	70	75	80	90	
HL	50	55	50	50	
EM	60	60	70	80	
Word List	NU6-50	NU6-25	NU6-50	NU6-25	
Method	CD	CD	Tape	CD	

LEFT EAR WORD RECOGNITION

	1	2	3	4	5
%	80				
HL	80				
EM	85				
Word List	NU6-25				
Method	CD				


Buttons: Save Changes, View Audiogram, Print, Help <?>, eXit Program

In the EFFECTIVE MASKING LEVELS section, INITIAL SRT refers to the first Speech Reception Threshold readings obtained during testing. Values obtained after masking, re-instructing the patient, or retesting for reliability purposes should be placed in the RETEST field. If data is in both fields, the retest data will be displayed on the graphic audiogram.

In the WORD RECOGNITION TESTING section, selections in the WORD LIST fields include both 25 and 50 word lists of the most commonly used lists.

If the cursor is paused over an editable field, a hint with acceptable values is displayed.

Whenever changes are saved, you will see a message similar to this one:



Acoustic Immittance Tab

Click on the Acoustic Immittance tab or enter ALT-I to make the Acoustic Immittance tab active.

This tab allows for entry of tympanometry and acoustic reflex information.

The screenshot shows a software window titled "Audiometric Data for ZZTESTPATIENT,PAMS(0140) as seen on JUN 18, 2003". It has four tabs: "Audiogram Entry", "Pure Tones", "Speech Audiometry", and "Acoustic Immittance", with the last one being active. The interface is divided into several sections for data entry:

- TYMPANOMETRY**: Contains fields for "Middle Ear Pressure" (range: -600 to 400 daPa), "Ear Canal Volume" (range: .1 to 6 cc), and "Peak Static Acoustic Immittance" (range: 0.01 to 10 mmhos). It has separate columns for "Probe Right" and "Probe Left".
- ACOUSTIC REFLEX THRESHOLDS**: This section is further divided into:
 - CONTRALATERAL**: Fields for "Stim Left" and "Stim Right" at 500 Hz, 1000 Hz, 2000 Hz, and 4000 Hz.
 - IPSILATERAL**: Similar fields for "Stim Right" and "Stim Left" at the same frequencies.
 - REFLEX DECAY**: Fields for "Stimulus Right" and "Stimulus Left" at 500 Hz and 1000 Hz, with a dropdown for "negative" or "positive".
 - REFLEX HALF-LIFE**: Fields for "Stimulus Right" and "Stimulus Left" at 500 Hz and 1000 Hz.

At the bottom of the window are five buttons: "Save Changes", "View Audiogram", "Print", "Help <?>", and "Exit Program".


If the cursor is paused over an editable field, a hint with acceptable values is displayed.

The Audiogram Display

1. Accessing the Audiogram Display

Typically, this application can be invoked from either CPRS or a stand-alone desktop shortcut (see instructions below). If either of these methods is desired but is not available, contact your facility's IRM Service to have the necessary installation or setup procedures completed.


When invoking the application, if your facility does not have single signon enabled, you may need to log in with your local **VistA** Access and Verify codes. A form similar to the following will appear.

A screenshot of a Windows-style dialog box titled "VISTA Sign-on". The main text area contains a message: "You are online with the Denver Distribution Center 'DEVELOPMENT' system. This system is available as a non-production simulated VistA environment for application development." The text is flanked by lines of asterisks. Below the text is a logo consisting of three interlocking 'V' shapes. To the right of the logo are two input fields: "Access Code:" followed by a field containing six asterisks, and "Verify Code:" followed by a field containing eight asterisks. To the right of these fields are two buttons: "OK" with a green checkmark icon and "Cancel" with a red 'X' icon. At the bottom of the dialog, there are four labels with corresponding text: "Server: AMMDDCDHC2", "Volume: TNR", "UCI: DEV", and "Port: //./nul:411".

In some cases, an end user may be set up for access to multiple broker environments. If that is the case, refer to [Appendix C](#) for further instructions.

Invoking From the Desktop:


Double-click the desktop icon for the *Audiogram Display* application. After logging in to your local broker system (if necessary), enter the beginning of the patients name for faster lookup.

A screenshot of a small dialog box titled "Find Patient...". It has a single text input field with the label "Enter start of last name:" above it. The text "ZZTESTPA" is entered into the field. Below the field are two buttons: "OK" and "Cancel".


Invoking From CPRS:

From the *Cover Sheet* in the CPRS application, click on the **T**ools menu. *Audiogram Display* should be one of the options on that menu. Since a patient has already been selected in CPRS, you will not have to choose one.


Then select from a list of audiograms in the Audiometric Exam Data file (#509850.9) that matches the name entered. (Be sure the selected line has the yellow dots before pressing OK.



If a match was not found, you will receive the following notice:




Depending on how the application was entered, you may see the list of entries in the AUDIOMETRIC EXAM file.



At this point you may select, or click on the Find button to narrow the search again.

If no record is available or you do not select to view one, you will receive the following alert message.



The application will then terminate.


This page intentionally left blank

2. Display Conventions

All data in the display examples in this document is test data only and may not reflect normal readings from real patients. See [Appendix A](#) for an explanation of how values from the record are selected to appear on the graph when two readings are taken.

All menu options at the top of the display are available with short-cut keys also. If the index character does not appear underlined automatically, underlining will appear by holding down the ALT key. Pressing the ALT key and the underlined letter key will trigger the option. They are also available by holding down the right mouse key and selecting the action.

Access to a Key, showing symbol meanings, is available from the option selections.



The Print All option will print the entire display, as viewable at the time the option is selected.

The frequency, in Hertz (Hz), is represented logarithmically on the horizontal axis (abscissa) in values from 125 to 8000 Hz. The dashed lines for 750, 1500, 3000, and 6000 Hz are placed on the graph in the logarithmic position.

The hearing level (HL), in decibels (dB), is represented on the vertical axis (ordinate) in values from -10 to 120. See [Appendix D](#) for further information.

Pure Tone symbols are drawn on the display so that the midpoint of the symbol centers on the vertical ruling and the horizontal axis at the appropriate hearing level.

Bone-conduction symbols are drawn adjacent to, but not touching the frequency coordinate ruling and centered vertically at the hearing level. The symbol for the left ear is placed to the right of the vertical ruling and that for the right ear is to the left of the vertical ruling.


By selecting the Copy I menu bar selection, the image referenced is copied to the Windows clipboard, from where it may be pasted into any picture-ready document (e.g., MS Word).

This page left intentionally blank

3. Main View

Audiogram in overlapping view

The default view shows the audiogram with overlapping right ear and left ear series. As is common practice, the right ear series displays in red and the left ear series in blue. Note the choices on the menu bar (see previous page on [Display Conventions](#)).




This page left intentionally blank

4. Separate View

The default view displays overlapping right- and left-ear measurements. Right and left readings may be separated for each audiogram displayed, as shown below. Use the *Separate I* menu bar selections to accomplish this, and the *Overlapping I* selections to return to the default view. All views display the date of the exam and age of the patient on that date above the audiogram. The name and title of the person who conducted the exam appear below the audiogram.

ZZTESTPATIENT,PAMS (609040140); On: Dec 23, 2002; Age: 62



This page left intentionally blank

5. Viewing 10-2364

The *Table 1* menu bar selection, presents a computer-generated VA Standard Form 10-2364 containing the values from the selected audiogram. The values in the form intended only for printing or viewing. The user may enter a short comment in the field provided. This comment is not saved to the database, but only for printing the current view.

VETERANS AFFAIRS AUDIOLOGICAL EVALUATION																			
AIR CONDUCTION																			
RIGHT Ear										LEFT Ear									
	250	500	1000	1500	2000	3000	4000	6000	8000		250	500	1000	1500	2000	3000	4000	6000	8000
Threshold	10	20	30*	35	40	45	50	45	60		30	25	35	40*	45	50	55	60	80
Masking																			
BONE CONDUCTION																			
Threshold			45*		45	50					55*+	65*		60	70				
Masking																			
ACOUSTIC IMMITTANCE										PURE TONE AVERAGES									
	Peak Pres	Vec	Static Meas			Peak Pres	Vec	Static Meas		Ear	2 freq	3 freq	4 freq						
Probe Right	200	2	3		Probe Left	300	3	4		RIGHT	25	30	41						
Contra AR Thresholds					Reflex Decay	Contra AR Thresholds					Reflex Decay	LEFT	30	35	46				
Pr-R	500	1000	2000	4000	500	1000	2000	4000	500	1000	Pr-L	500	1000	2000	4000	500	1000		
St-L	70	80	90	100	Neg	Pos	St-R	50	60+	75	95	Neg	Pos						
Ipsi AR Thresholds					Half-Life	Ipsi AR Thresholds					Half-Life	INTER-TEST-CONSISTENCY							
St-R	75				5	6	St-L	80				7	8	RIGHT	LEFT				
										FAIR			GOOD						
SPEECH AUDIOMETRY																			
SRT			RIGHT Ear Word Recognition										LEFT Ear Word Recognition						
Ear	1	2		1	2	3	4	5	PBMAX	PI/PB		1	2	3	4	5	PBMAX	PI/PB	
R	50	40	%	30	35	40			40	0.13	%	60					60		
L	60	80	HL	90	80	75					HL	80							
ML	60	90	ML		80						ML	80							
Comments:																			
Patient Name [Last,First MI] ZZTESTPATIENT,PAMS										Age 62		SSN 609040140							
Examining Audiologist URRUTIA,PAM										Date of Exam Dec 23, 2002									
Examining Station or Clinic DENVER (DDC)																			
Computer generated VA form 10-2364 AUDIOLOGICAL EVALUATION Developed at Denver Distribution Center for QUASAR 3.0*3 2003																			

The options to *copy*, *print* or *exit* this form are available by clicking on the *right mouse button*, while the cursor is over the form.

This page left intentionally blank

Glossary

Acronyms

AC	Air Conduction
AR	Acoustic Reflex
ASPS	Audiology & Speech Pathology Service
BC	Bone Conduction
C&P	Compensation and Pension
CAR	Contralateral Acoustic Reflex
CNC	Consonant Nucleus Consonant
dB	Decibel
DDC	Denver Distribution Center
EM	Effective Masking
HL	Hearing Level
Hz	Hertz
IAR	Ipsilateral Acoustic Reflex
MCL	Most Comfortable Loudness
ML	Masking Level
NR	No Response
NU	Northwestern University
PI/PB	Performance Intensity-Phonetically Balanced words
Pr-L	Probe Left
Pr-R	Probe Right
PSAS	Prosthetics & Sensory Aids Service
PTA	Pure Tone Average (see Appendix D)
SAT	Speech Awareness Threshold
SRT	Speech Reception Threshold
UCL	Uncomfortable Loudness

This page left intentionally blank

Appendices

A: Determination of Series Values for Display

Some of the key rules applied in preparing series values for display in the *Audiogram Display* application are as follows:

Both the initial and final thresholds and their respective masking levels are obtained.

If the repeat threshold indicated a 'No Response' (+), there will be a gap in the series.

The initial threshold value for the series will be replaced with the repeat value for any of the following reasons:

If the repeat value is masked and the initial is not.

If the repeat value is not masked, and the repeat is less than the initial.

If the initial reading does not have a value.

B: Calculation Of PB MAX And PI/PB

Calculation of PB MAX and PI/PB values in the *Audiometric Exam Enter/Edit* application is based on specific industry-standard formulae established for these measurements. Basic descriptions of some of these formulae are as follows:

PB MAX is the maximum percentage from the word recognition testing.

PB MIN is the minimum percentage from the word recognition testing.


PI/PB is an indice of possible retrocochlear pathology. It is calculated from the formula: $(PB\ MAX - PB\ MIN) / (PB\ MAX)$.

The PI/PB index will assess multiple scores and levels in one ear. The calculation will occur only when a second score obtained at a higher presentation level, is poorer than a prior score at a lower presentation level.

The result is always a value less than 1.0.

C: If You Have Access To Multiple Broker Environments

If you have additional access to a training or development broker environment, you may need to select the environment you want to use from a form like the following:



By selecting the correct server name from the drop down list, you will access the correct environment. (Using the down arrow key will also move from one selection to another)

See your local IRM or ADPAC to obtain the correct server name.

D: Calculation of Pure Tone Averages

Pure Tone Averages displayed in the *Audiometric Exam Enter/Edit* application are automatically calculated and supplied when sufficient information is entered. At each level, the numbers used are the ones that would appear on the graph for that level.

The formulas used are:

PTA TWO = Average of the lowest two readings from 500, 1000, and 2000 Hz.

PTA THREE = Average of 500, 1000 and 2000 Hz.

PTA FOUR = Average of 1000, 2000, 3000 and 4000 Hz.

E: VA FileMan Date/time Formats

Examples of Valid Dates:

JAN 20 1957 or 20 JAN 57 or 1/20/57 or 012057

T (for TODAY), T+1 (for TOMORROW), T+2, T+7, etc.

T-1 (for YESTERDAY), T-3W (for 3 WEEKS AGO), etc.

If the year is omitted, the computer uses CURRENT YEAR.

Two-digit year assumes no more than 20 years in the future, or 80 years in the past.

If only the time is entered, the current date is assumed.

Follow the date with a time, such as JAN 20@10, T@10AM, 10:30, etc.

You may enter a time, such as NOON, MIDNIGHT or NOW.

You may enter NOW+3' (for current date and time Plus 3 minutes

*Note--the Apostrophe following the number of minutes)

Note: Time is REQUIRED in this response.

This page left intentionally blank

Index

10-2364, 25

Acoustic Immittance, 14

AGE AT VISIT, 8

Audience, i

Benefits, ii

Bone-conduction symbols, 19

Calculation Of PB MAX And PI/PB, 30

central database, ii

changes, 13

Conventions, 1, 2, 19

Date Signed, 11

DATE/TIME OF VISIT, 8

default, 21

delete, 9

Determining Series Values, 29

disabled fields, 2

EXAMINING AUDIOLOGIST, 8

formula, 30

frequency, 19

hearing level, 19

horizontal axis, 19

Hospital Location, 8

Key, 19

logarithmic scale, 19

Masking level fields, 2, 10

menu options, 19

PB MAX, 30

PB MIN, 30

PI/PB, 30

Print All, 19

Pure Tone Averages, 11, 12

Pure Tone symbols, 19

Purpose, i

REFERRAL SOURCE, 8

ROES, ii

select the system, 30

short-cut key combinations, 1

single signon, 4, 15

Speech Audiometry, 12

SRT, 12

symbol meanings, 19

tab stops, 2

Table, 25

TYPE OF VISIT, 9

VA ELIGIBILITY STATUS, 8, 9

vertical axis, 19

View, 3

word lists, 13