Chapter 7: Viewing and Printing BCMA Reports

Viewing and Printing BCMA Reports	.7-3
Benefits of This Chapter	.7-3
Enhancements to Report Printing.	.7-3
PRN Effectiveness List	
MAH Date Range Selection	7-3
Future Orders	7-3
Viewing/Printing Demographics Data	.7-4
Viewing/Printing Allergy and ADR Information	.7-6
Viewing/Printing an Order from Inpatient Medications	.7-8
Viewing/Printing a Due List Report	
Viewing/Printing a Medication Log Report	.7-16
Viewing/Printing an MAH Report	7-19
Viewing/Printing a Missed Medications Report	7-22
Viewing/Printing a PRN Effectiveness List Report	.7-25
Viewing/Printing an Administration Times Report	
Viewing/Printing a Medication History Report	

Benefits of This Chapter

This chapter describes a variety of reports originating from BCMA that you can now print (and yes, still view on-screen) via the VDL or via Tabs in CPRS. You can also use it to view patient demographic and allergy information.

Enhancements to Report Printing

You now have the option of printing hard copy reports by Patient or by Ward — directly from the VDL — or to continue viewing the reports on-screen as was the case in Version 1.0. When you select a Ward, BCMA lets you choose whether to print a Ward Report by Patient or by Room-Bed. BCMA automatically provides you with a list of printers available for outputting your reports.



TIP:

All reports from BCMA print in a 132-column format. Here's a few other printing enhancements included with this new version of BCMA:

PRN Effectiveness List

You can now print the PRN Effectiveness List from the Reports menu. In the previous version, you could only print this report using the CHUI version of BCMA.

MAH Date Range Selection

Now you can specify the date range (as defined by CPRS) that you want to view or print the MAH Report. If no parameter is defined, the range defaults to a seven-day range, as in the previous version of BCMA.

Future Orders

Another exciting feature is the ability to view (or print) future active medication orders using the Due List command in the Reports menu. You can determine whether to include active Unit Dose and/or IV medication orders on the Report, plus Future Orders and related Addendums/Changes.

Viewing/Printing Demographics Data



TIP:


Use the Patient Demographics Button in the upper left-hand corner of the VDL, below the Tool Bar, to quickly display patient demographics data. You can use the Patient Demographics command (or button) to view or print personal, admission, eligibility, and appointment information about the patient whose orders are displayed on the VDL. This information was electronically documented when the patient was admitted to your medical center.

To View/Print a Patient's Demographics Data:

1 Select the Patient Demographics command from the View menu. The Patient Inquiry dialog box displays, with the patient's demographics data.

Keyboard Only Users: Press **ALT+V** to display the View menu, and then press **P** to display the Patient Inquiry dialog box.

Example: Patient Inquiry Dialog Box



Viewing/Printing Demographics Data (cont.)

To View/Print a Patient's Demographics Data: (cont.)

- **2** Perform one of the following actions:
 - Review the patient's demographics data, and then click CANCEL to return to the patient's VDL.
 - Click **PRINT** to display the Printer dialog box. Proceed to step #3.

Keyboard Only Users: Use **TAB** to move among the **PRINT** and **CANCEL** buttons.

3 Select a printer from the drop-down list box that you want to use for outputting a report with the patient's demographics data. Then click **OK**. An Information message displays.

Note: The printer that you select becomes the "default" printer for *all* reports printed from BCMA.

- 4 Note the task number for your print job, and then click **OK** to return to the Patient Inquiry dialog box. At the dialog box, click **CANCEL** to return to the patient's VDL.
- **5** Retrieve the hard copy report, containing the patient's demographics data, from your printer.

Viewing/Printing Allergy and ADR Information



TIP:

Use the Allergies
Button on the
Tool Bar to quickly
display an
alphabetical listing
of all food, drug,
and "other" allergy
and ADR
information
documented about
the patient.

Clicking on the Allergies command (or button) displays allergy and adverse drug reaction information documented about the patient in the ART package. This includes the causative agent, drug class, signs and symptoms, and whether the allergy/ADR was verified and observed. If the Allergies command/button is grayed out, no allergies or adverse drug reactions were documented about the patient.


Note: BCMA only displays allergy information, not a drug "interaction" check.

To View/Print a Patient's Allergy and ADR Information:

1 Select the Allergies command from View menu. The Patient Allergy List dialog box displays, with the patient's allergy and ADR information.

Keyboard Only Users: Press **ALT+V** to display the View menu, and then press **A** to display the Patient Allergy List dialog box.

Example: Patient Allergy List Dialog Box



Viewing/Printing Allergy and ADR Information (cont.)

To View/Print a Patient's Allergy and ADR Information: (cont.)

- **2** Perform one of the following actions:
 - Review the patient's allergy and ADR information, and then click **CANCEL** to return to the patient's VDL.
 - Click **PRINT** to display the Printer dialog box. Proceed to step #3.

Keyboard Only Users: Use **TAB** to move among the **PRINT** and **CANCEL** buttons.

3 Select a printer from the drop-down list box that you want to use for outputting a report with the patient's allergy and ADR information. Then click **OK**. An Information message displays.

Note: The printer that you select becomes the "default" printer for *all* reports printed from BCMA.

- 4 Note the task number for your print job, and then click **OK** to return to the Patient Allergy List dialog box. At the dialog box, click **CANCEL** to return to the patient's VDL.
- **5** Retrieve the hard copy report, containing the patient's allergy and ADR information, from your printer.

Viewing/Printing an Order from Inpatient Medications



TIP:

Double-click on a medication order, displayed on the VDL, to view the details of the order from Inpatient Medications.


BCMA now provides several ways for you to view (and now print) Unit Dose and IV medication orders from Inpatient Medications V. 5.0. They include the following: 1) highlight a medication on the VDL, and then press **F4**; 2) double-click on an order displayed on the VDL; 3) use the Display Order command in the Due List menu; or 4) use the Display Order command in the Right Click drop-down menu.

To View/Print an Order from Inpatient Medications:

- 1 Highlight a medication order on the patient's VDL.
- **2** Select the Display Order command from the Due List menu. The Display Order dialog box displays, with the details of the medication order.

Keyboard Only Users: Press **F4** to display the Display Order dialog box.


Example: Display Order Dialog Box for Unit Dose Medication Order




Viewing/Printing an Order from Inpatient Medications (cont.)

This section provides examples of an IV Piggyback and IV medication order as they were entered using Inpatient Medications V. 5.0.

Example: Display Order Dialog Box for IV Piggyback Medication Order



Example: Display Order Dialog Box for IV Medication Order



Viewing/Printing an Order from Inpatient Medications (cont.)

To View/Print an Order from Inpatient Medications: (cont.)

- **3** Perform one of the following actions:
 - Review the patient's medication order, and then click **CANCEL** to return to the patient's VDL.
 - Click **PRINT** to display the Printer dialog box. Proceed to step #4.

Keyboard Only Users: Use **TAB** to move among the **PRINT** and **CANCEL** buttons.

4 Select a printer from the drop-down list box that you want to use for outputting a report, with the patient's medication order. Then click **OK**. An Information message displays.

Note: The printer that you select becomes the "default" printer for *all* reports printed from BCMA.

- 5 Note the task number for your print job, and then click **OK** to return to the Display Order dialog box. At the dialog box, click **CANCEL** to return to the patient's VDL.
- **6** Retrieve the hard copy report, containing the patient's medication order information, from your printer.

Viewing/Printing a Due List Report



TIP:

BCMA lists information on the Due List Report by IV orders, then by Unit Dose orders. Within each section, medications are listed alphabetically.

BCMA now lets you print (and yes, still view) a Due List Report directly from the VDL. This report provides detailed information about active *and* future Unit Dose and IV medication orders that are "due" for administering to a patient — during a timeframe that you specify — within a 24-hour period.


The Due List Report includes patient demographic data, allergy and ADR information, plus detailed information about an order, such as whether (or *not*) the medication is a self-med; the medication type, schedule, dose, and route; Special Instructions; administration times; Last Given date and time; Start/Stop date and time; and the individual(s) who verified the order.

To View/Print a Due List Report:

1 Select the Due List command from the Reports menu. The Patient Due List dialog box displays.

Keyboard Only Users: Press ALT+R to display the Reports menu, and then press **D** to display the Patient Due List dialog box.

Example: Patient Due List Dialog Box





TIP:

The Patient Due
List dialog box
'defaults' include all
check boxes
selected, and
today's date listed.

Viewing/Printing a Due List Report (cont.)



TIP:

You can type the information in the list boxes, or use the drop-down arrows to make your selection.



TIP:

The bottom of the Due List Report includes blank lines for making any "Changes/
Addendums" to the patient's orders.
Future orders display on a separate page.

To View/Print a Due List Report: (cont.)

2 Use the **DOWN ARROW**, within the list boxes, to select the Date, and Start and Stop Times of the Due List Report that you want to view on-screen or print.

Note: Clicking on the **DOWN ARROW** on the Date for Report list box displays a calendar. You can use the scroll arrows in the upper corners of the calendar to display a different month, and then click on a date to select it and close the calendar. Clicking inside the Start/Stop Time list boxes provides selections.

Keyboard Only Users: Use **TAB** to move among the different areas of the dialog box.

3 In the Schedule Types, Order Types, and Include areas, click inside a check box to include your selection on the Due List Report.

Keyboard Only Users: Press **SPACEBAR** to select (check) a check box.

4 In the Print by area, click inside a Radio button to print the Due List Report by Patient or by Ward.

Note: If you choose to print the Due List Report by Ward/Room-Bed, make your selection from the list box provided. This Report lists information alphabetically by patient.

Keyboard Only Users: Use the **DOWN ARROW** to select the Ward Radio button, and a ward location from the drop-down list box.

- **5** Click **OK** to display the Due List Report on-screen.
- **6** Perform one of the following actions:
 - Review the Due List Report, and then click **CANCEL** to return to the patient's VDL.
 - Click **PRINT** to display the Printer dialog box. Proceed to step #7.

Keyboard Only Users: Use **TAB** to move among the **PRINT** and **CANCEL** buttons.

7 Select a printer from the drop-down list box that you want to use for outputting the Due List Report. Then click **OK**. An Information message displays.


Note: The printer that you select becomes the "default" printer for *all* reports printed from BCMA.

Viewing/Printing a Due List Report (cont.)

To View/Print a Due List Report: (cont.)

- 8 Note the task number for your print job, and then click **OK** to return to the Patient Due List dialog box. At the dialog box, click **CANCEL** to return to the patient's VDL.
- **9** Retrieve the Due List Report from your printer. An example of each report type is provided in this section.


Example: Due List Report By Patient



Viewing/Printing a Due List Report (cont.)

This section provides an example of a Due List Report by Ward/Patient.


Example: Due List Report By Ward/Patient



Viewing/Printing a Due List Report (cont.)

This section provides an example of a Due List Report by Ward/Room-Bed.

Example: Due List Report By Ward/Room-Bed



Viewing/Printing a Medication Log Report



TIP:

Medications ordered and documented using the CPRS Med Order Button in BCMA display on the Med Log Report opposite the text titled "BCMA/CPRS Interface Entry."



TIP:

Use the Med Log Button on the BCMA Tool Bar, or the Reports Tab in CPRS, to quickly display the Patient Medication Log dialog box. You can now print (and still view on-screen) a patient's Medication Log Report (or Med Log Report), listing *every* action taken on a medication order within a specified 24-hour period. You can choose to include Comments and Audits performed on the patient's medication orders. Audits include "actions" taken on a medication.

The Medication Log Report also includes patient demographic data, allergy and ADR information, plus detailed information about the order, such as the drug/additive/solution; the orderable item dose, schedule, route, and injection site; Action date and time; Activity Start/Stop date and time; Unique Identifier Number for the bag; Unit ordered and given, and the clinician(s) who took an action on the medication.


Note: The Medication Log Report lists medications with a "Not Given" status in the Audit Trail section of the Log, *not* on the VDL. Administrations that were Cancelled are *not* listed on the Report.

To View/Print a Medication Log Report:

1 Select the Medication Log command from the Reports menu. The Patient Medication Log dialog box displays.

Keyboard Only Users: Press **ALT+R** to display the Reports menu, and then press **L** to display the Patient Medication Log dialog box.





Viewing/Printing a Medication Log Report (cont.)



TIP:

You can type the information in the list boxes, or use the drop-down arrows to make your selection.

To View/Print a Medication Log Report: (cont.)

2 Use the **DOWN ARROW**, within the list boxes, to select the Date, and Start and Stop Times of the Medication Log Report that you want to view on-screen or print.

Note: Clicking on the **DOWN ARROW** on the Date for Report list box displays a calendar. You can use the scroll arrows in the upper corners of the calendar to display a different month, and then click on a date to select it and close the calendar. Clicking inside the Start/Stop Time list boxes provides selections.

Keyboard Only Users: Use **TAB** to move among the different areas of the dialog box.

3 In the Include area, click inside a check box to include your selection on the Medication Log Report.

Keyboard Only Users: Press **SPACEBAR** to select (check) a check box.

4 In the Print by area, click inside a Radio button to print the Medication Log Report by Patient or by Ward.

Note: If you choose to print the Medication Log Report by Ward/Room-Bed, make your selection from the list box provided. This Report lists information alphabetically by patient.

Keyboard Only Users: Use the **DOWN ARROW** to select the Ward Radio button, and a ward location from the drop-down list box.

- 5 Click **OK** to display the Medication Log Report on-screen.
- **6** Perform one of the following actions:
 - Review the Medication Log Report, and then click **CANCEL** to return to the patient's VDL.
 - Click **PRINT** to display the Printer dialog box. Proceed to step #7.

Keyboard Only Users: Use **TAB** to move among the **PRINT** and **CANCEL** buttons.

7 Select a printer from the drop-down list box that you want to use for outputting the Medication Log Report. Then click **OK**. An Information message displays.

Note: The printer that you select becomes the "default" printer for *all* reports printed from BCMA.



TIP:


See the section
"Viewing/Printing a
Due List Report" to
see how a Ward
Report by RoomBed looks when
printed.

Viewing/Printing a Medication Log Report (cont.)

To View/Print a Medication Log Report: (cont.)

- **8** Note the task number for your print job, and then click **OK** to return to the Patient Medication Log dialog box. At the dialog box, click **CANCEL** to return to the patient's VDL.
- **9** Retrieve the Medication Log Report from your printer. An example is provided below.

Example: Medication Log Report By Patient



Unique

Identifier

Number for IV Bag

Viewing/Printing an MAH Report



TIP:

You can specify the date range (as defined for CPRS), for an MAH Report. See "Note" on this page for more information.

TIP:



Use the Med Admin
Hist Button on the
BCMA Tool Bar, or
the Reports Tab in
CPRS, to quickly
display the Patient
Medications Given
dialog box.

You can now view on-screen *and* print an MAH Report for Unit Dose and IV medication orders. This Report lists a clinician's name and initials, and the exact time that an action was taken on an order (in a conventional MAR format). Each order is listed alphabetically by the orderable item. The date column lists three asterisks (***) if a medication was Discontinued.

An MAH Report includes patient demographic data, allergy and ADR information, plus detailed information about the order, such as the drug/additive/solution; the medication schedule, dose, route, and injection site; the actual Administration Times; the name and initials of the clinician who administered the medication; and the individuals who verified the order. It also includes information about when an order is placed on, and taken off Hold by a Provider.


Note: If no parameter is defined in CPRS, the maximum date range defaults to a seven-date range, as in the previous version of BCMA. For example, a Report would list the Sunday proceeding, and the Saturday following, the date that you selected for the Report.

To View/Print an MAH Report:

1 Select the Medication Admin History command from the Reports menu. The Patient Medications Given dialog box displays.

Keyboard Only Users: Press **ALT+R** to display the Reports menu, and then press **H** to display the Patient Medications Given dialog box.

Example: Patient Medications Given Dialog Box



Viewing/Printing an MAH Report (cont.)



TIP:

You can type the information in the list boxes, or use the drop-down arrows to make your selection.

To View/Print an MAH Report: (cont.)

2 Use the **DOWN ARROW**, within the list boxes, to select the Start and Stop Dates of the MAH Report that you want to view on-screen or print.

Note: Clicking on the **DOWN ARROW** on the Start and Stop Date list boxes displays a calendar. You can use the scroll arrows in the upper corners of the calendar to display a different month, and then click on a date to select it and close the calendar. Clicking inside the Start/Stop Time list boxes provides selections.

Keyboard Only Users: Use **TAB** to move among the different areas of the dialog box.

3 In the Print by area, click inside a Radio button to print the MAH Report by Patient or by Ward.

Note: If you choose to print the MAH Report by Ward/Room-Bed, make your selection from the list box provided. This Report lists information alphabetically by patient.

Keyboard Only Users: Use the **DOWN ARROW** to select the Ward Radio button, and a ward location from the drop-down list box.

- 4 Click **OK** to display the MAH Report on-screen.
- **5** Perform one of the following actions:
 - Review the MAH Report, and then click **CANCEL** to return to the patient's VDL.
 - Click **PRINT** to display the Printer dialog box. Proceed to step #6.

Keyboard Only Users: Use **TAB** to move among the **PRINT** and **CANCEL** buttons.

6 Select a printer from the drop-down list box that you want to use for outputting the MAH Report. Then click **OK**. An Information message displays.

Note: The printer that you select becomes the "default" printer for *all* reports printed from BCMA.

7 Note the task number for your print job, and then click **OK** to return to the Patient Medications Given dialog box. At the dialog box, click **CANCEL** to return to the patient's VDL.



TIP:

See the section
"Viewing/Printing a
Due List Report" to
see how a Ward
Report by RoomBed looks when
printed.

Viewing/Printing an MAH Report (cont.)



TIP:


You can print this report, after a patient has been discharged, using CHUI BCMA.

To View/Print an MAH Report: (cont.)

8 Retrieve the MAH Report from your printer. An example is provided below.

Note: If you select dates that do *not* have applicable data, the MAH Report will appear "empty" when displayed on-screen or printed.

Example: MAH Report By Patient



Viewing/Printing a Missed Medications Report



TIP:

Medications placed on "Hold" in CPRS or Inpatient Medications V. 5.0, appear on the Missed Medications Report with the word "Hold" in parentheses to the right of them. BCMA V. 2.0 lets you now print a Missed Medications Report (or Missed Meds Report), for Continuous and One-Time Unit Dose and IV Piggyback medications that were *not* administered to a patient during a Med Pass.

A Missed Medications Report includes patient demographic data, allergy and ADR information, plus detailed information about the order, such as the medication type; the administration date and time; and the order number.


Note: Self-medications no longer display on the Missed Medications Report.

To View/Print a Missed Medications Report:

1 Select the Missed Medications command from the Reports menu. The Patient Missed Medications dialog box displays.

Keyboard Only Users: Press **ALT+R** to display the Reports menu, and then press **M** to display the Patient Missed Medications dialog box.

Example: Patient Missed Medications Dialog Box



Viewing/Printing a Missed Medications Report (cont.)



TIP:

You can type the information in the list boxes, or use the drop-down arrows to make your selection.



See the section
"Viewing/Printing a
Due List Report" to
see how a Ward
Report by RoomBed looks when
printed.

To View/Print a Missed Medications Report: (cont.)

2 Use the **DOWN ARROW**, within the list boxes, to select the Date, and Start and Stop Times for the Missed Medications Report that you want to view on-screen or print.

Note: Clicking on the **DOWN ARROW** on the Date for Report list box displays a calendar. You can use the scroll arrows in the upper corners of the calendar to display a different month, and then click on a date to select it and close the calendar. Clicking inside the Start/Stop Time list boxes provide selections.

Keyboard Only Users: Use **TAB** to move among the different areas of the dialog box.

3 In the Print by area, click inside a Radio button to print the Missed Medications Report by Patient or by Ward.

Note: If you choose to print the Report by Ward/Room-Bed, make your selection from the list box provided. This Report lists information alphabetically by patient.

Keyboard Only Users: Use the **DOWN ARROW** to select the Ward Radio button, and a ward location from the drop-down list box.

- 4 Click **OK** to display the Missed Medications Report on-screen.
- **5** Perform one of the following actions:
 - Review the Missed Medications Report, and then click **CANCEL** to return to the patient's VDL.
 - Click **PRINT** to display the Printer dialog box. Proceed to step #6.

Keyboard Only Users: Use **TAB** to move among the **PRINT** and **CANCEL** buttons.

6 Select a printer from the drop-down list box that you want to use for outputting the Missed Medications Report. Then click **OK**. An Information message displays.

Note: The printer that you select becomes the "default" printer for *all* reports printed from BCMA.

7 Note the task number for your print job, and then click **OK** to return to the Patient Missed Medications dialog box. At the dialog box, click **CANCEL** to return to the patient's VDL.


Viewing/Printing a Missed Medications Report (cont.)

To View/Print a Missed Medications Report: (cont.)

8 Retrieve the Missed Medications Report from your printer. An example is provided below.

Note: The column "Order Num" on the Missed Medications Report lists the actual order number of the medication from Inpatient Medications V. 5.0. This information is quite helpful when troubleshooting problems with BCMA.

Example: Missed Medications Report By Patient



Viewing/Printing a PRN Effectiveness List Report



TIP:


You can print this report, after a patient has been discharged, using CHUI BCMA. Now you can view (and print) a PRN Effectiveness List Report using BCMA V. 2.0. This Report lists PRN medications administered to a patient that need Effectiveness comments. It also includes patient demographic data, allergy and ADR information; plus the PRN medication, administration date and time; and the individual(s) who administered the order.

To View/Print a PRN Effectiveness Report:

1 Select the PRN Effectiveness List command from the Reports menu. The PRN Effectiveness List dialog box displays.

Keyboard Only Users: Press **ALT+R** to display the Reports menu, and then press **P** to display the PRN Effectiveness List dialog box.

Example: PRN Effectiveness List Dialog Box



Viewing/Printing a PRN Effectiveness List Report (cont.)



TIP:

You can type the information in the list boxes, or use the drop-down arrows to make your selection.

To View/Print a PRN Effectiveness List Report: (cont.)

2 Use the **DOWN ARROW**, within the list boxes, to select the Date, and Start and Stop Times for the PRN Effectiveness List Report that you want to view on-screen or print.

Note: Clicking on the **DOWN ARROW** on the Date for Report list box displays a calendar. You can use the scroll arrows in the upper corners of the calendar to display a different month, and then click on a date to select it and close the calendar. Clicking inside the Start/Stop Time list boxes provide selections.

Keyboard Only Users: Use **TAB** to move among the different areas of the dialog box.

3 In the Print by area, click inside a Radio button to print the PRN Effectiveness List Report by Patient or by Ward.

Note: If you choose to print the Report by Ward/Room-Bed, make your selection from the list box provided. This Report lists information alphabetically by patient.

Keyboard Only Users: Use the **DOWN ARROW** to select the Ward Radio button, and a ward location from the drop-down list box.

- 4 Click **OK** to display the PRN Effectiveness List Report on-screen.
- **5** Perform one of the following actions:
 - Review the PRN Effectiveness List Report, and then click **CANCEL** to return to the patient's VDL.
 - Click **PRINT** to display the Printer dialog box. Proceed to step #6.

Keyboard Only Users: Use **TAB** to move among the **PRINT** and **CANCEL** buttons.

6 Select a printer from the drop-down list box that you want to use for outputting the PRN Effectiveness List Report. Then click **OK**. An Information message displays.

Note: The printer that you select becomes the "default" printer for *all* reports printed from BCMA.

7 Note the task number for your print job, and then click **OK** to return to the PRN Effectiveness List dialog box. At the dialog box, click **CANCEL** to return to the patient's VDL.



TIP:


See the section
"Viewing/Printing a
Due List Report" to
see how a Ward
Report by RoomBed looks when
printed.

Viewing/Printing a PRN Effectiveness List Report (cont.)

To View/Print a PRN Effectiveness List Report: (cont.)

8 Retrieve the PRN Effectiveness List Report from your printer. An example is provided below.

Example: PRN Effectiveness List Report By Patient



Viewing/Printing an Administration Times Report

You can now print, *and* still view on-screen, an Administration Times Report. This report lists a patient's medications by the scheduled administration time (from the earliest to the latest).

The Administration Times Report includes patient demographic data, allergy and ADR information, plus detailed information about the order, such as the medication type, dose, and route; and the administration time.


Note: Self-medications no longer display on the Administration Times Report.

To View/Print an Administration Times Report:

1 Select the Administration Times command from the Reports menu. The Patient Ward Administration Times dialog box displays.

Keyboard Only Users: Press **ALT+R** to display the Reports menu, and then press **A** to display the Patient Ward Administration Times dialog box.

Example: Patient Ward Administration Times Dialog Box



Viewing/Printing an Administration Times Report (cont.)



TIP:

You can type the information in the list boxes, or use the drop-down arrows to make your selection.

To View/Print an Administration Times Report: (cont.)

2 Use the **DOWN ARROW**, within the list box, to select the date of the Administration Times Report that you want to view on-screen or print.

Note: Clicking on the **DOWN ARROW** on the Date for Report list box displays a calendar. You can use the scroll arrows in the upper corners of the calendar to display a different month, and then click on a date to select it and close the calendar.

Keyboard Only Users: Use **TAB** to move among the different areas of the dialog box.

3 In the Print by area, click inside a Radio button to print the Administration Times Report by Patient or by Ward.

Note: If you choose to print the Administration Times Report by Ward/Room-Bed, make your selection from the list box provided. This Report lists information chronologically by patient.

Keyboard Only Users: Use the **DOWN ARROW** to select the Ward Radio button, and a ward location from the drop-down list box.

- 4 Click **OK** to display the Administration Times Report on-screen
- **5** Perform one of the following actions:
 - Review the Administration Times Report, and then click **CANCEL** to return to the patient's VDL.
 - Click **PRINT** to display the Printer dialog box. Proceed to step #6.

Keyboard Only Users: Use **TAB** to move among the **PRINT** and **CANCEL** buttons.

6 Select a printer from the drop-down list box that you want to use for outputting the Administration Times Report. Then click **OK**. An Information message displays.

Note: The printer that you select becomes the "default" printer for *all* reports printed from BCMA.

7 Note the task number for your print job, and then click **OK** to return to the Patient Ward Administration Times dialog box. At the dialog box, click **CANCEL** to return to the patient's VDL.



TIP:


See the section
"Viewing/Printing a
Due List Report" to
see how a Ward
Report by RoomBed looks when
printed.

Viewing/Printing an Administration Times Report (cont.)

To View/Print an Administration Times Report: (cont.)

8 Retrieve the Administration Times Report from your printer. An example is provided below.

Example: Administration Times Report By Patient



Viewing/Printing a Medication History Report



TIP:

Right click on a medication displayed on the VDL, and then select the Med History command in the Right Click drop-down menu to display the Medication History dialog box.

BCMA now provides two ways to access a patient's Medication History Report (or Med History Report). They include the Med History command in the Due List menu *and* in the Right Click drop-down menu.


This Report lists the administration date and time, and orderable item of a medication highlighted on the VDL. It also includes the medication status, schedule type, and dose; room location; and initials of the clinician who administered the medication.

Note: A Medication History Report is called an "Administration History Report" in CPRS. You can access it by selecting the Meds Tab in CPRS, and then right clicking on a medication.

To View/Print a Medication History Report:

1 Select the Med History command from the Due List menu. The Medication History dialog box displays.

Keyboard Only Users: Press **ALT+D** to display the Due List menu, and then press **E** to display the Medication History dialog box.



Example: Medication History Dialog Box

Viewing/Printing a Medication History Report (cont.)

To View/Print a Medication History Report: (cont.)

- **2** Perform one of the following actions:
 - Review the patient's Medication History Report, and then click **CANCEL** to return to the patient's VDL.
 - Click **PRINT** to display the Printer dialog box. Proceed to step #3.

Keyboard Only Users: Use **TAB** to move among the **PRINT** and **CANCEL** buttons.

3 Select a printer from the drop-down list box that you want to use for outputting the Medication History Report for the medication highlighted on the VDL. Then click **OK**. An Information message displays.

Note: The printer that you select becomes the "default" printer for *all* reports printed from BCMA.

- 4 Note the task number for your print job, and then click **OK** to return to the Medication History dialog box. At the dialog box, click **CANCEL** to return to the patient's VDL.
- **5** Retrieve the Medication History Report from your printer.