

# **SURGERY**

## **TECHNICAL MANUAL**

Version 3.0 July 1993

(Revised April 2003)

## **Revision History**

The table below lists changes made since the initial release of this manual. Use the Change Pages document to update an existing manual or use the entire updated manual.

**Note:** The Change Pages document may include unedited pages needed for two-sided copying. Only edited pages display the patch number and revision date in the page footer.

Date	Revised Pages	Patch Number	Description
04/03	Title Page, i, v – vi, 67-69	SR*3*115	<ul> <li>- Updated Title page, Revision History, and Table of Contents.</li> <li>- Added new fields TIME OUT VERIFIED and IMAGING CONFIRMED to the SURGERY file (#30) listing; this addition caused the pages to wrap.</li> </ul>
02/03	Title Page, I, 15, 69-70, 71, 124, 185	SR*3*107	<ul> <li>- Updated Title page and Revision History.</li> <li>- Added new field AUTOMATED CASE CART ORDERING to the SURGERY SITE PARAMETERS file (#133) listing.</li> <li>- Added new fields SPD COMMENTS and DYNAMED NOTIFIED to the SURGERY file (#130) listing; this addition caused the pages to wrap.</li> <li>- Added callable routine CSLSUR1 to the Calls Made by Surgery section.</li> </ul>
12/02	185 - 186	SR*3*111	- Added callable routine DGUTL4 to the Calls Made by Surgery section.
11/02	Title Page, i, (ii) 13, (14) (55), 56 (119), 120 (185), 186	SR*3*109	Updated Revision Date and Revision History.  - Changed the field name DEFAULT BLOOD REQUEST to DEFAULT BLOOD COMPONENT.  - Changed the definition of the field REQ BLOOD KIND from a pointer to the BLOOD PRODUCT file (#66) to a free text field.  - Deleted DEFAULT BLOOD REQUEST field, and replaced with DEFAULT BLOOD COMPONENT field.  - Added callable routine VBECA5A to the Calls Made by Surgery section.
07/93			Original Release of Technical Manual.

<This page intentionally left blank for two-sided printing.>

# **Table of Contents**

Introduction	
Special Notations	
Change Pages	
Package Management	
How This Software is Organized	
-	
Menu Outline	
List Scheduled Operations	6
Implementation and Maintenance	13
Installation	
Site Specific Parameters	13
Security Keys	23
File Security	27
•	
Files	
SURGERY (130)	
PERSON FIELD RESTRICTION (131)	
SURGERY TRANSPORTATION DEVICES (131.01)	
OPERATION TIMES (131.25)	
SURGERY DISPOSITION (131.6)	
OPERATING ROOM (131.7)	
SURGERY UTILIZATION (131.8)	
PROSTHESIS (131.9)	
SURGERY POSITION (132)RESTRAINTS AND POSITIONAL AIDS (132.05)	110
SURGICAL DELAY (132.4)	
ANESTHESIA SUPERVISOR CODES (132.95)	117
SURGERY SITE PARAMETERS (133)	
SURGERY INTERFACE PARAMETER (133.2)	
MONITORS (133.4)	
IRRIGATION (133.6)	
SURGERY REPLACEMENT FLUIDS (133.7)	
SURGERY WAITING LIST (133.8)	
OPERATING ROOM TYPE (134)	
SURGERY CANCELLATION REASON (135)	120
SKIN PREP AGENTS (135.1)	
SKIN INTEGRITY (135.2)	
PATIENT MOOD (135.3)	
PATIENT CONSCIOUSNESS (135.4)	
SURGERY TRANSCRIPTION (136)	
PERIOPERATIVE OCCURRENCE CATEGORY (136.5)	
(100.0)	

## Contents

LOCAL SURGICAL SPECIALTY (137.45)	132
ELECTROGROUND POSITIONS (138)	133
RISK MODEL LAB TEST (139.2)	133
Options	138
Option Descriptions	138
Exported Menu Structure	161
Routines	163
Callable Routines	165
Cross References	167
Archiving and Purging	183
Archiving	183
Purging	183
External Relations	185
Packages Needed to Run Surgery	185
Calls Made by Surgery	
Data Base Integration Agreements (DBIAs)	187
Internal Relations	188
Generating On-line Documentation	190
Additional Information	192
Templates	
Input Templates	192
Sort Template	192
Editing Input Templates	193
Globals in the Surgery Package	193
Journaling Globals	
Routine Mapping	
Resetting DTIME	194
Package-Wide Variables	196
Glossary	200

### 63 IV STARTED BY

**POINTER** 

This is the name of the person that started the IV for this operative procedure.

64 CULTURES

WORD-PROCESSING

These are the names of cultures sent to the laboratory for examination.

### 65 SURGERY POSITION

POINTER

This is the position in which the patient is placed for this operative procedure. This information will appear on the Nurse Intraoperative Report.

.01 SURGERY POSITION

**POINTER** 

This is the position in which the patient is placed for this operative procedure. More than one position may be entered for each case.

1 TIME PLACED

DATE/TIME

This is the date/time that the patient was placed in this position. Times without a date can be entered.

66 PRIN DIAGNOSIS CODE

**POINTER** 

This is the principal ICD9 diagnosis code. It should be entered for all cases and will be used for Surgery Central Office reporting needs.

67 CANCELLATION AVOIDABLE

SET

This field contains a set of codes used to flag a cancellation as being avoidable or unavoidable. It is used when determining the percentage of avoidable cancellations.

68 SCHEDULED PROCEDURE

FREE TEXT

This field contains the scheduled (or original) principal procedure for this case. It will be compared to the actual procedure completed.

69 CODING VERIFIER

POINTER

This is the person who last updated procedure and/or diagnosis descriptions and/or codes for this case using the Update/Verify Procedure/Diagnosis Codes [SRCODING EDIT] option. This field is updated automatically by the option when information is changed.

70 CANCELLED BY

POINTER

This is the name of the person who cancelled this surgical case. This information is automatically entered when a case is cancelled.

71 TIME OUT VERIFIED

SET

This field refers to the completion of a "Time Out" verification process prior to the start of the procedure. Enter YES if the "Time Out" verification process was completed prior to the start of the procedure. If entered "NO", a justification should be documented in the Nursing Care Comments.

72 IMAGING CONFIRMED

SET

This field refers to the completion of the verification process for the presence of relevant imaging data to confirm that the operative site for the correct patient are available, properly labeled and properly presented, and verified by two members of the operating team prior to the start of the procedure.

### 80 SPD COMMENTS

WORD-PROCESSING

This field contains any information for SPD that cannot be entered elsewhere. These comments will be sent to SPD via the Surgery/CoreFLS interface.

### 81 DYNAMED NOTIFIED

SET

This field indicates whether or not a notification has been sent to DynaMed by way of the CoreFLS interface. YES indicates at least one notification has been sent, while a null value or zero indicates that no notification has been sent. The first notification sent to DynaMed will be a NEW APPOINTMENT notification. Subsequent notifications will be for editing, canceling, or deleting notifications, as appropriate.

### 100 ORDER NUMBER

**POINTER** 

This is the pointer to the ORDER file (100). It will be created when a case is requested.

### 101 STAFF/RESIDENT

SET

This determines whether the surgeon was a resident or staff. It will be used for categorizing procedures in the Annual Report of Surgical Procedures.

### 102 REASON FOR NO ASSESSMENT

SET

This is the reason why no assessment was entered for this particular surgical case. It should be entered if any major procedure was excluded from the risk assessment module.

- 1 Patient did NOT receive general, spinal, or epidural anesthesia.
- 2 Number of surgical cases entered into the Surgical Risk Study exceeded 36 over an 8-day time frame.
- 3 Number of TURPs or TURBTs exceeded 5 cases over an 8-day time frame.
- 4 Study exclusion criteria prohibits patient entry.
- 6 Surgical Clinical Nurse Reviewer was on Annual Leave.
- 8 Case was a concurrent case, secondary to an assessed primary case.
- 9 Number of inguinal hernias exceeded 5 cases over an 8-day time frame.

### 103 ANESTHETIST CATEGORY

SET

This field holds the category of the principal anesthetist, which is used on the Anesthesia AMIS report to enumerate the number of anesthetics administered by each category.

### 118 NON-OR PROCEDURE

SET

This field is a flag signifying this case is a non-OR surgical procedure.

### 119 NON-OR LOCATION

POINTER

This is the location (file 44) where this non-OR procedure was performed.

### 120 DATE OF PROCEDURE

DATE/TIME

This is the date that the non-OR procedure was performed. The date of procedure must be entered for all non-OR cases.

- TIME PROCEDURE BEGAN DATE/TIME
  This is the date and time that the non-OR procedure began.
- TIME PROCEDURE ENDED DATE/TIME
  This is the date and time that all the procedures for this non-OR case are complete.
- PROVIDER

  This is the person who performs the major portion of the principal non-OR procedure. This field is required for several reports.
- 124 ATTEND PROVIDER POINTER
  This is the name of the attending staff provider responsible for this case.
  This information appears on several reports.
- 125 MEDICAL SPECIALTY POINTER
  This is the medical specialty credited for doing this non-OR procedure. Many reports are sorted by the medical specialty. This field should be entered prior to completion of this non-OR procedure.
- PROCEDURE OCCURRENCE FREE TEXT
  This is an occurrence that is related to a non-O.R. procedure. If there are not any non-O.R. procedure occurrences, this field should be left blank. Do not enter 'NO' or 'NONE'.
  - .01 PROCEDURE OCCURRENCE FREE TEXT This is an occurrence that is related to a non-O.R. procedure. If there are not any non-O.R. procedure occurrences, this field should be left blank. Do not enter 'NO' or 'NONE'.
  - OUTCOME TO DATE SET
    This is the outcome to date of this non-O.R. procedure occurrence.
  - 2 DATE OCCURRENCE NOTED DATE/TIME
    This is the date that this occurrence was noted. The time of day can be entered, but is not required.
  - TREATMENT INSTITUTED FREE TEXT
    This is the type of treatment instituted as a result of this non-O.R. procedure occurrence.
  - 4 OCCURRENCE COMMENTS WORD-PROCESSING This is information that may be helpful in documentation of the non-O.R. procedure occurrence.
  - 5 OCCURRENCE CATEGORY POINTER
    This is the name of the category for which this occurrence will be grouped for Surgery Central Office reporting needs.

# 127 SEQUENTIAL COMPRESSION DEVICE SET This determines whether a sequential compression device was used.

### 128 LASER TYPE

### FREE TEXT

This determines whether a laser was used during this procedure. If applicable, enter the type of laser used during this surgical procedure.

### 200 OPERATIONS THIS ADMISSION

### **NUMERIC**

This is the total number of surgical procedures, prior to the index or principal operation, which required the patient to be taken to the operating room for any type of surgical intervention during this hospital admission. Include all procedures whether or not they are part of the inclusion/exclusion criteria.

### 201 REDO PROCEDURE

#### SET

This determines whether the principal operative procedure was a reoperation in the same anatomic location for the same purpose as the first operation regardless of the length of time from the original surgical date.

### 202 CURRENT SMOKER

#### SET

For non-cardiac assessment, did the patient smoke cigarettes in the year prior to admission for surgery?

For cardiac assessment, did the patient smoke tobacco in any form in the 2 weeks prior to surgery?

### 202.1 PACK/YEARS

### NUMERIC

If the patient has ever been a smoker, enter the total number of pack/years of smoking for this patient. The number of pack/years is determined by multiplying the number of packs of cigarettes smoked per day by the number of years the patient has smoked. If the patient has never been a smoker, enter "0". If the smoking history for this patient cannot be determined, enter "NS". The possible range for number of pack/years is 0 to 200.

### 203 HISTORY OF COPD

### SET

This determines whether the patient has chronic obstructive pulmonary disease (COPD) resulting in functional disability, and/or hospitalization, and/or requiring chronic bronchodilator therapy, and/or an FEV1 of less than 75% of predicted.

For non-cardiac assessment, do not include patients with acute asthma, an acute and chronic inflammatory disease of the airways resulting in bronchospasm.

### 204 VENTILATOR DEPENDENT

SET